


LIQUIBAND[®] XL
SKIN CLOSURE SYSTEM

www.liquiband.com

PROTECTION

WITHOUT COMPROMISE




 **LIQUIBAND[®]**
TOPICAL ADHESIVES

PROTECTION

WITHOUT COMPROMISE

Designed to meet both patient and clinician needs

LiquiBand® XL is a product combining both a mesh patch and a liquid adhesive to securely hold large and high-tension incisions closed for up to 14 days.*



Signature winged applicator

- ▶ Allows simple, easy, one-handed priming and adhesive application⁷

Self-adhering hexagonal mesh

- ▶ Easy to apply and holds skin edges closed while applying liquid adhesive

Water-resistant film

- ▶ Allows patients to shower with physician approval^{1,*}

Microbial barrier as long as the film remains intact^{2,3}

- ▶ Provides protection from Gram-positive, Gram-negative, and fungal microbes

Consistent dry time in as little as 60 seconds⁴

- ▶ Faster than the leading market competitor⁶

Non-invasive skin closure system^{4,5}

- ▶ Does not puncture skin like sutures or staples do, providing cosmesis without track marks

2-octyl cyanoacrylate adhesive formulation

- ▶ Durable and effective for up to 14 days




PROTECTION

WITHOUT COMPROMISE

Strength you expect...

LiquiBand® XL provides equivalent wound holding strength to other large wound closure systems^{8, 9} and is stronger than 3-0 USP sutures^{3, 9, 10}


Wound Closure Strength




...in less time

LiquiBand® XL dries faster than the leading large wound closure systems⁶ and securely closes incisions faster than sutures.¹¹

Drying Times



Closure Times



LIQUIBAND® XL

SKIN CLOSURE SYSTEM

PROTECTION WITHOUT COMPROMISE



DESCRIPTION	QTY PER BOX	PRODUCT CODE
22cm Mesh (1 x 22cm x 4cm Mesh Strip , 1 x 2.7g Adhesive Applicator)	2	72014019
44cm Mesh (2 x 22cm x 4cm Mesh Strips, 2 x 2.7g Adhesive Applicators)	2	72014027
66cm Mesh (3 x 22cm x 4cm Mesh Strips, 3 x 2.7g Adhesive Applicators)	2	72014046

Advanced Medical Solutions Limited

Western Wood Way, Langleigh Science Park, Plymouth, Devon, PL7 5BG, United Kingdom

📍 www.liquiband.com | ✉ info@liquiband.com | ☎ +44 (0) 8444 125 754



Advanced Medical Solutions

References: Data on File at Advanced Medical Solutions Ltd.

* Please refer to the IFU for complete application instructions.

© Based on *in vitro* testing of 3-0 and 4-0 USP polyamide sutures

LiquiBand® is a registered trademark of Advanced Medical Solutions Limited. All Rights Reserved.

Copyright © Advanced Medical Solutions Limited 2023

IRM 06 0325 03